

ASSODONT

Carrello CR3A



Carrello CR3A con 3 ripiani arrotondati

Carrello elettrificato, in acciaio elettrozincato, verniciato con polveri epoxidiche. Si compone di 3 eleganti ripiani arrotondati con bordo in tubolare per ostacolare la caduta degli oggetti.

Ciò che rende il carrello veramente unico è la sua funzionalità e l'accurato design. Maneggevolezza, spigoli arrotondati, raffinatezza delle finiture, dimensioni ridotte, si accompagnano alla praticità delle 2 prese elettriche che fanno del carrello non solo una pratica base di appoggio, ma un utilissimo strumento per lo spostamento e l'utilizzo di piccole apparecchiature.

CR3A trolley with 3-rounded-shelves

Electrified trolley with 3 wealthy rounded-shelves with tubular edges, 4 casters and 2 useful general purpose sockets for electrical appliances.

It is made of epoxy powder painted electro-zink plated steel and tubular supports.

What makes the trolley unique is its functionality and the accuracy of its design. Handiness, round edges, accurate finishes and small dimensions are the main characteristics of this trolley. It is also the ideal and rational solution to the small problems that may arise in the organization of a dental surgery.

Colore standard	RAL 9002	Basic colour
Ingombro	Mm L=400 x P=450 x H=800	Max overall dimensions
Dimensioni ripiani	Mm L= 400 x P= 300	Shelves dimensions
Portata max per ripiano	Kg 10	Max weight capability per shelf
Peso	Kg 15	Weight
Imballo	Mm L=500 x P=360 x H=760	Packaging